BILANCIA DIGITALE PESA BAMBINI DIGITAL BABY SCALE

INMD-039
MANUALE D'USO
USER MANUAL

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Thank you for purchasing the INMD-039 baby weighing scale.

Before using the product, it is important to carefully read the instructions. We recommend keeping this manual for future reference.

WARNINGS AND IMPORTANT SAFETY INSTRUCTIONS

- Before installing and using the device, carefully read the instruction manual and keep it for future reference. Incorrect or improper use of the device may damage the scale and cause harm to the user.
- Never leave the baby unattended during the weighing process.
- Use the scale only for weights below 20 kg. Using the scale with higher weights may damage the device and be dangerous.
- To clean the scale, use a soft cloth, possibly dampened with water. Do not immerse it in liquids or rinse it under running water.
- Do not use aggressive detergents/solvents or abrasive substances for cleaning.
- While using the scale, keep it away from sources of electromagnetic interference such as microwave ovens, mobile phones, and wireless network devices.
- Do not short-circuit the battery contacts inside the battery compartment.
- In case of malfunction, do not attempt to repair the scale

- yourself. Always contact qualified technical personnel.
- Never allow the baby to play with the packaging materials, especially plastic bags, which may cause suffocation.
 - If a malfunction occurs due to electrostatic charges, remove the batteries for a few minutes. This will restore the normal functioning of the device.
 - Do not place anything on the scale when it is not in use.
 - Do not expose the scale to prolonged loads, as this may damage the measuring mechanism.
 - Removing or replacing the batteries will delete all recorded weight measurements.
 - The scale will automatically turn off after 25 seconds of inactivity. To turn off the scale immediately, press and hold the "ON/OFF/TARE" button until the scale turns off.
 - Avoid impacts or falls and place the scale on a flat, stable, and secure surface.
 - This scale is intended for home use only and is designed to measure the baby's weight. It is not approved for medical use or any purpose other than private use. Since the scale is exclusively for private use, it is exempt from legal periodic accuracy checks required for commercial transactions.

TECHNICAL SPECIFICATIONS

- Large surface: 50.6 x 23.3 cm
- Product dimensions: 54.8 x 32 x 3.5 cm
- LCD display: 58 x 25 mm
- Capacity: 20 kg
- Graduation: 10 g
- Units of measurement: kg/lb/st
- Tare function
- Automatic shut-off
- Overload indication: O-Ld
- Low battery indication: Lo
- Power supply: 2 AAA batteries (not included)
- Product weight: 1.17 kg

OPERATION

1. Open the battery compartment under the scale and insert 2 AAA batteries.
2. Place the scale on a hard, flat surface.
3. Press the button **ON/OFF
TARE** to turn on the scale.
4. Do not load or weigh anything until the display resets to zero. The scale will be ready for use when the LCD screen shows "0.00 kg."
5. Press the UNIT button to select the desired unit of measurement.
6. Place the baby on the weighing platform; the display will show the weight.
7. Tare function: Press the button **ON/OFF
TARE**.
8. The maximum capacity is 20 kg. If "O-Ld" appears on the display, it means the weight exceeds the scale's capacity.
9. If the low battery indication "Lo" appears, replace the batteries promptly. Ensure they are inserted with the correct polarity.
10. If no further weighing is needed, press the button **ON/OFF
TARE** or wait for the automatic shutdown in a few seconds.

NOTES

1. Do not immerse the scale in water or clean it with corrosive chemicals.
2. When weighing, place the scale on a hard, flat surface. Do not use it on carpets or soft materials.

RATING LABEL

Innoliving Spa
 Via Merloni 2/B - 60131 Ancona - Italy
 INMD-039 BILANCIA DIGITALE PESA BAMBINI
 DIGITAL BABY SCALE
 MAX: 20 Kg DIV: 10g 3V == (2xAAA batt.)
 SOLO USO DOMESTICO - HOUSEHOLD USE ONLY

LOT n. 
 MADE IN CHINA 

This product complies with all applicable European directives.

INFORMATION TO USERS according to Legislative Decree No. 49 of March 14, 2014 " 2012/19/UE Implementation of the Directive on Waste Electrical and Electronic Equipment (WEEE) " The crossed-out dustbin symbol indicates that the product at the end of its life must be

 collected separately from other waste. The user should, therefore, take the equipment with the essential components at the end of its useful life to the separate collection center of electronic and electrical waste, or return it to the retailer when purchasing new equipment of equivalent type, in ratio of one to one, or one to zero for devices with larger side less than 25 CM. The separate collection for the decommissioned equipment for recycling, treatment and environmentally compatible disposal contributes to avoid possible negative effects on the environment and human health and promotes recycling of the materials. Improper disposal of the product by the user entails the application of administrative sanctions according to Legislative Decree No. 49 of 14 March 2014.

The batteries used in this device must be disposed of in the special bins at the end of their life. Please inform yourself about the local rules on separate collection of batteries. The correct disposal of batteries helps preventing potentially negative consequences on the environment and human health.

Innoliving Spa
Via Merloni, 2/B
60131 Ancona Italy
Tel 071 2133550
www.innoliving.it

MADE IN CHINA

Rev.00_02.2025